


# Q.E.D.<sup>®</sup>

SALIVA ALCOHOL TEST



*The On-The-Spot Leader In  
Alcohol Testing . . .*



- Simple
- Accurate
- Reliable
  - D.O.T. Approved
  - CLIA Waived
  - *In Vitro* Diagnostic Use



The disposable Q.E.D.<sup>®</sup> Saliva Alcohol Test is a revolutionary breakthrough in quantitative on-site alcohol detection.

Distributed by:

# Q.E.D.<sup>®</sup>

SALIVA ALCOHOL TEST

1-800-869-3538

## Three easy steps:



1. Swab mouth to collect saliva.



2. Insert collector into test.



3. Read color bar.




**OraSure Technologies, Inc.**

diagnostic solutions for the new millennium

220 East First Street • Bethlehem, PA 18015  
www.orasure.com

## On-Site Simplicity with Clinical Lab Accuracy . . .



The Q.E.D.<sup>®</sup> test can be used as an accurate screen in place of blood tests. In clinical trials, saliva alcohol levels measured with the Q.E.D.<sup>®</sup> tests demonstrated high correlation to blood analyzed by gas chromatography ( $r=0.98$ )<sup>1</sup> and produced no false positive or false negative results<sup>2</sup>.

Data on file. <sup>1</sup>Study conducted by Thomas Jefferson University Hospital, Department of Emergency Medicine.  
<sup>2</sup>Study conducted by Rutgers University, Center of Alcohol Studies.

## It's Reliable, Fast, Simple . . .

- Built-in QA Spot™ assures quality control for every test
- Quantitative results in minutes
- Long shelf life
- Reads like a thermometer
- No instrumentation
- Non-invasive

## The Q.E.D.<sup>®</sup> test is used for . . .

- Pre-employment or pre-duty
- Random testing
- Post-accident
- Reasonable cause
- Follow-up testing
- Testing upon admission
- Interventions
- Counseling sessions
- Random testing
- Pre-discharge
- Admission screening
- Routine trauma treatment protocol
- Seizure or unconscious victims

**Worksite testing  
for DOT or drug-free  
workplace programs**

**Alcohol treatment  
centers**

**Emergency/trauma  
departments**

## Additional Information . . .

- Available testing range:  
→A150 (0-145 mg/dL, 0-0.145%)
- DOT Screening Test Technician Certification available from Nationwide Network of Trainers or **NEW** DOT-approved training video
- Toll Free Technical Support—Monday thru Friday 8 a.m. - 6 p.m. (eastern time)

### OraSure Technologies

Product Number	Description
31150(B,C)	Q.E.D. <sup>®</sup> A150—case of 30 tests(B), case of 100 tests(C)
31050S	Q.E.D. <sup>®</sup> DOT Ethanol Control (also 31150S, 31350S)
31150T	Screening Test Technician (STT) Training Video Kit
31150X	Additional STT Training Student Workbooks
31000V	Q.E.D. <sup>®</sup> Procedure Video—case of 10